Vakuumtechnik im System
Technology for Vacuum Systems
Technologie du vide

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© 2001 VACUUBRAND GMBH + CO KG Printed in Germany

D-97866 Wertheim - Postfach 1664
D-97877 Wertheim - Alfred-Zippe-Str. 4
Tel. (+49)9342 / 808-0 - Fax (+49)9342 / 59880

Contents

- ➔ **Safety information!**
- ➔ **Technical data**
- ➔ **Use and operation**
- ➔ **Troubleshooting**
- ➔ **Readjustment of the vacuum gauge**
- ➔ **Calibration in the factory**
- ➔ **Replacement of the battery**
- ➔ **Notes on return to the factory**
Repair - Return - DKD-Calibration
- ➔ **Health and safety clearance form**

Safety information!

Remove all packing material, remove the product from its packing-box, remove the protective covers and keep them, inspect the equipment.

If the equipment is damaged, notify the supplier and the carrier in writing within three days; state the item number of the product together with the order number and the supplier's invoice number. Retain all packing material for inspection.

Do not use the equipment if it is damaged.

If the equipment is not used immediately, replace the protective covers. Store the equipment in suitable conditions.

- ☞ Read this manual before installing or operating the equipment.
- ☞ Obey all applicable regulations and realize all required measures.
- ☞ Position the vacuum gauge and the vacuum connection in the vacuum system in a way that no draining of condensate towards the pressure transducer inside the vacuum gauge can occur.

The current supply of the devices is realized by the incorporated 9V battery.

- ☞ Recommended type see section "Technical data".
- ☞ The battery is free from mercury and cadmium.
- ☞ The battery must not be charged, opened, thrown into fire, heated to more than 100°C or brought in contact with water.
- ☞ When replacing the battery, always insert it correctly into the battery compartment.
- ☞ When using a battery with a different chemical system (e. g. alkaline), lifetime may be considerably shorter.

- ☞ The use of an accumulator is not recommended because of its low capacity.

The devices are **not suitable** when working with **dangerous or explosive gases or explosive or flammable mixtures**. Ensure that the materials of the wetted parts are compatible (see "Technical data"). If necessary adopt suitable measures.

Max. permitted pressure at the vacuum gauge: 2 bar absolute.

- ☞ The display flashes at pressures above 1060 mbar.

Max. ambient temperature for storage and operation: 60°C.

Max. temperature of gaseous media at the pressure transducer (gas!): 80°C.

Ensure that maintenance is done only by suitably trained and supervised technicians. Ensure that the maintenance technician is familiar with the safety procedures which relate to the equipment processed by the vacuum system and that the equipment if necessary is approximately decontaminated before starting maintenance.

Use only **genuine spare parts and accessories**.

- ☞ Otherwise safety and performance of the equipment as well as the electromagnetic compatibility of the equipment might be reduced.

In order to comply with law (occupational, health and safety regulations and regulations for environmental protection) products returned to the manufacturer can be **repaired / DKD calibrated** only when following certain procedures (see section "Notes on return to the factory").

Technical data

| Type | DVR 2 |
|--|---|
| Pressure transducer | capacitive / ceramic |
| Automatic sleep mode | user selectable operation time 1-999 minutes or "On" (continuous operation) (factory-set: 15 minutes) |
| Pressure reading | analogue und digital LCD display mbar, Torr or hPa (to be switched between) |
| Measuring cycle time | selectable: automatic, 1 x within 3s, 1 x within 1s, 3 x within 1s (factory-set: automatic) |
| Measuring range | 1 mbar - 1080 mbar (1 Torr - 810 Torr) |
| Max. permitted pressure at the vacuum gauge | 2 bar absolute |
| Uncertainty (with transducer carefully adjusted and at constant temperature) | <+/-1 mbar (0.75 Torr) +/-1 digit |
| Temperature coefficient | <+/- 0.07 mbar/K (<0.05 Torr/K) |
| Max. temperature of gaseous media at the pressure transducer (gas!) | +80°C |
| Permitted ambient temperature range (operation) | +10°C to +60°C |
| Permitted temperature range (storage) | -10°C to +60°C |
| Vacuum connection | small flange NW 16 and screw-in stepped hose nozzle for vacuum hoses with 6/10 mm inside diameter |
| Power supply | battery 9V lithium, 1.2 Ah, Ultralife U 9VL |
| Mass | 375 g |
| Dimensions of housing (L x W x H) | 115 mm x 115 mm x 66 mm |

| Components | Wetted parts |
|-----------------------------|------------------------|
| Vacuum connection | PTB |
| Seals | FPM (e. g. Viton®*) |
| Pressure transducer housing | stainless steel |
| Pressure transducer | aluminiumoxide-ceramic |

* reg. trade mark Du Pont

We reserve the right for technical modification without prior notice!


Use and operation

Display


adjustment mode/
warning symbol

automatic sleep mode


replace battery/
display of battery capacity

vacuum units
(to be switched between)

Keys


key ON/OFF

key UP/DOWN


key MODE

The keys are on the **rear side** of the device underneath the symbols at the front side.


Switch between vacuum units:

- Switch off. Press key MODE simultaneously with key ON/OFF.
- ☞ The vacuum unit as from the last operation is displayed.


- Select desired vacuum unit by using key UP/DOWN.

- Press key ON/OFF to confirm the selected vacuum unit and to terminate mode.


Pressure measurement:

Upon switching on, the absolute pressure (total pressure, independent of gas type) ranging from 1 mbar to 1080 mbar is displayed in the selected vacuum unit (Torr, mbar or hPa).

- ☞ The pressure transducer has been adjusted in the factory. In general, there is no need for readjustment by the user because of the excellent longtime stability (see section "Readjustment").


Clock symbol


- ☞ The vacuum gauge is equipped with an automatic sleep mode to save battery power when the instrument is not used. When the clock symbol appears on the display, the vacuum gauge will switch off automatically after approx. 1 minute.
- ➔ Press key ON/OFF, the clock symbol disappears and the preselected operation time will start again.

Adjusting the operating time and measuring cycle

- ☞ The time until **automatic switching off (operating time)** can be set by the user within a range of 1 to 999 minutes (factory-set: 15 minutes) or **continuous operation "On"**.


- ➔ Press key MODE to select the time setting mode.
- ☞ The clock symbol and the actually set operating time in minutes are displayed.


- ☞ The arrow in the right direction indicates that the time can be increased by using key UP/DOWN.
- ➔ Press key UP/DOWN to increase operating time.

To decrease operating time:


- ➔ Press key MODE to change arrow direction to the left.
- ➔ Press key UP/DOWN to decrease operating time.
- ➔ Press key ON/OFF to confirm.
- ☞ The device switches to mode "**adjustment of measuring cycle**".
- ☞ The **measuring cycle** can be chosen between Automatic, 1 measurement in 3s, 1 measurement in 1s and 3 measurements in 1s.

To adjust measuring cycle:


- ☞ The measuring cycle as from the last operation (factory-set: "A") is displayed.
- ➔ Press key UP/DOWN to select the measuring cycle: A = Automatic (automatic determination of the measuring cycle; frequent measurement in case of great pressure change), S1 = 1 measurement in 3s, S2 = 1 measurement in 1s, S3 = 3 measurements in 1s.
- ☞ The flashing of the arrow (only "S1", "S2" or "S3") symbolizes the frequency of the pressure measurement).
- ☞ Press key MODE to change the direction of the adjustment.

- ➔ Press key ON/OFF to confirm the measuring cycle and to terminate mode.

Display of the status of discharge of the battery

When switching off the device, the remaining battery capacity is displayed qualitatively by the number of arrows in the display for a few seconds.

- ☞ If the battery symbol appears in the display during measurement, the battery should be replaced.
- ☞ The correct function of the device is no longer ensured.

Troubleshooting

| Fault | Possible cause | Remedy |
|--|---|---|
| <input type="checkbox"/> No display or display disappears. | <ul style="list-style-type: none"> ➔ Battery empty? ➔ Automatic sleep mode? ➔ Other cause? | <ul style="list-style-type: none"> ☞ Replace battery. ☞ Switch on, increase operating time if necessary. ☞ Contact local distributor. |
| <input type="checkbox"/> Incorrect display. | <ul style="list-style-type: none"> ➔ Pressure transducer has drifted off? ➔ Moisture in the pressure transducer? ➔ Battery empty, battery symbol is displayed? | <ul style="list-style-type: none"> ☞ Readjust pressure transducer. ☞ Dry, respectively evacuate the transducer and readjust if necessary. ☞ Replace battery. |
| <input type="checkbox"/> Display is flickering and/or battery symbol appears. | <ul style="list-style-type: none"> ➔ Battery empty? | <ul style="list-style-type: none"> ☞ Replace battery. |
| <input type="checkbox"/> Adjustment mode cannot be activated. | <ul style="list-style-type: none"> ➔ At the current pressure, readjustment is not possible. | <ul style="list-style-type: none"> ☞ Admit air to atmospheric pressure or evacuate to a pressure < 20 mbar. |
| <input type="checkbox"/> All segments of the LCD are displayed or no display although battery has been replaced. | <ul style="list-style-type: none"> ➔ Pressure transducer or measuring electronics defective? | <ul style="list-style-type: none"> ☞ Contact local distributor. |

Readjustment of the vacuum gauge


The equipment has been adjusted at the factory. In general, there is no need for adjustment by the user because of the excellent longtime stability.

Depending on operation conditions, type of application and accuracy requirements, an inspection and readjustment may become necessary.

Adjustment at atmospheric pressure

Admit air to the vacuum gauge. Make sure that the vacuum connection at the vacuum gauge is at atmospheric pressure.

Attention: Determine the exact actual atmospheric pressure, e. g. by using an accurate barometer or get accurate reading from the weather service, e. g. at the next airport, etc. (take into account the difference in altitude between e. g. airport and laboratory).


- Press key UP/DOWN simultaneously with key ON/OFF.

- ☞ The vacuum gauge switches to the adjustment mode (indicated by a warning triangle).

- ☞ The arrow in the right direction indicates that the pressure value can be increased by using key UP/DOWN.


- Press key UP/DOWN to increase the reading to actual atmospheric pressure.

To reduce the reading:

- Press key MODE to change arrow direction to the left.


- Then press key UP/DOWN to reduce the reading.

- Press key ON/OFF to confirm adjustment and to terminate mode.

Adjustment under vacuum

Evacuate the vacuum gauge to a pressure <0.5 mbar (<0.4 Torr) e. g. by applying a good rotary vane pump.


- Press key UP/DOWN simultaneously with key ON/OFF.

- ☞ The vacuum gauge switches to the adjustment mode (indicated by a warning triangle)

- ☞ The reading is adjusted automatically to "zero".


Attention: Adjustment under vacuum with a pressure higher than 0.5 mbar (0.4 Torr) reduces the accuracy of measurement. If the pressure is significantly higher than 0.5 mbar (0.4 Torr), adjustment to a reference pressure is recommended.

- Press key ON/OFF to confirm adjustment and to terminate mode.

Adjustment at reference pressure

Evacuate the vacuum gauge to an exact known reference pressure within the range of 0 20 mbar (0 15 Torr). Switch vacuum gauge to the adjustment mode (see "Adjustment under vacuum").

- Press key UP/DOWN to adjust the display from "0" ("zero") to the actual reference pressure in the vacuum line in the range of 0 to 20 mbar (0 to 15 Torr).
- Press key ON/OFF to confirm adjustment and to terminate mode.

Attention: The accuracy of the value of the reference pressure will directly affect the accuracy of the adjustment. If the nominal ultimate vacuum of a diaphragm pump is used as "reference" vacuum, the accuracy of the adjustment of the vacuum gauge might be doubtful. The diaphragm pump may not achieve the specified value due to condensate, poor state, failure of the valves or the diaphragm.

Calibration in the factory

Control of measuring equipment

The **VACUUBRAND DKD calibration laboratory** is accredited by the Physikalisch-Technische Bundesanstalt (PTB; German national institute for science and technology and the highest technical authority of the Federal Republic of Germany for the field of meteorology and certain sectors of safety engineering) for the measurable variable **pressure in the pressure range from 10^{-3} mbar to 1000 mbar** in accordance with the general criteria for the operation of testing laboratories defined in the standard EN ISO/IEC 17025:2000. Calibration in the VACUUBRAND calibration laboratory:

- To meet the requirements of the DIN ISO 9000ff and 10012 series of standards regarding the calibration of inspection, measuring and test equipment at specified intervals.
- To document that the vacuum gauges calibrated are traceable to the national standard of the PTB.

DKD recalibration **90 02 17**


Replacement of the battery

Note: Some of the components on the printed circuit board are sensitive to electrostatic discharge. Do not touch components. If necessary, ground the printed circuit board and the person replacing the battery suitably.

Attention: Use only a battery of the recommended type (see "Technical data"). The battery must not be charged or connected to an external voltage.

- Unscrew the four countersunk head screws at the rear side of the housing by using a Phillips screw driver.
- Disassemble the housing parts carefully.

- Remove the battery.
- Insert the new battery in the battery compartment correctly matching the positive and negative polarity (as shown by the engraving inside the compartment).
- ⚠ When reassembling the vacuum gauge, make sure that the switch pad and the printed circuit board are positioned correctly.
- Carefully screw both parts of the housing together.

Notes on return to the factory

Repair - return - DKD calibration

If the measures listed in section "Troubleshooting" do not work: Return the equipment to the factory for inspection and repair.

Opening or disassembly of the equipment voids all warranties and the licence if necessary.

Safety and health of our staff, laws and regulations regarding the handling of dangerous goods, occupational health and safety regulations and regulations regarding safe disposal of waste require that for all pumps and other products the "**Health and safety clearance form**" must be sent to our office duly completed and signed before any equipment is dispatched to our premises.

Fax or post a completed copy of the health and safety clearance form to us in advance. The declaration must arrive before the equipment. Enclose a second completed copy with the product. If the equipment is contaminated you must notify the carrier.

No repair / DKD calibration is possible unless the correctly completed form is returned. Inevitably, there will be a delay in processing the equipment if information is missing or if this procedure is not obeyed. Eventually, acceptance of the returned equipment will be refused.

If the product has come in contact with chemicals, radioactive substances or other substances dangerous to health or environment, the product must be decontaminated **prior to sending it back to the factory.**

To expedite repair and to reduce costs, please enclose a detailed description of the problem and the products' operation conditions with every product returned for repair. We submit **quotations** only on request and always at the customer's expense. If an order is given, the costs incurred are offset from the costs for repair or from the purchase price, if the customer prefers to buy a new product instead of repairing the defective one.

☛ **If you do not wish a repair on the basis of our quotation, the equipment might be returned to you disassembled and at your charge!**

We replace parts due to optical aspects upon your request (at your expense).

Before returning the equipment ensure that (if applicable):

- ☛ Equipment has been cleaned and/or decontaminated.
- ☛ All ports have been sealed.
- ☛ Equipment has been properly packed, if necessary, please order an original packaging (costs will be charged), marked as appropriate and the carrier has been notified.
- ☛ Ensure that the completed health and safety declaration is enclosed.

We hope for your understanding for these measures, which are beyond our control.

Scrapping and waste disposal:

Dispose of the equipment and any components removed from it safely in accordance with all local and national safety and environmental requirements. Particular care must be taken with components which have been contaminated with dangerous substances from the process. Do not incinerate fluoroelastomer seals and "O"-rings.

☛ You may authorize us to dispose of the equipment **at your expense.**

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Disclaimer: Our technical literature is only intended to inform our customer. The validity of general empirical values and results obtained under test conditions for specific applications depend on a number of factors beyond our control. It is therefore strictly the users' responsibility to very carefully check the validity of application to their specific requirements. No claims arising from the information provided in this literature will, consequently, be entertained.

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© 2002 VACUUBRAND GMBH + CO KG Printed in Germany 99 90 01

D-97866 Wertheim - Postfach 1664
D-97877 Wertheim - Alfred-Zippe-Str. 4
Tel. (++49)9342 / 808-0 - Fax (++49)9342 / 59880