

LIMITED OFFER!

SHOW
DEMO
UNIT



**Web Only
SALE!**
**While Supplies
Last**

To Order
Call 1-800-843-1794
or Email

customer-service@cglifesciences.com

Use promo code:
CG-14042



Where Chemistry Meets Biology

Microplate Genies

Units have aggressive mixing speeds from 750 to 3200 RPM for complete and uniform mixing in any microplate format. Small vortexing orbit of 1.0mm for thorough mixing regardless of sample viscosity. Speed and orbit combine to offer true vortexing action in each well of the microplate. Units accept most microplate types within the recommendations of the Society for Biomolecular Screening (SBS), even in 384-well formats. The Microplate Genie may be used with any single microplate. May be used in cold rooms and incubators for temperature sensitive applications. The built-in timer operates from 1 to 60 minutes or continuous.

Specifications:

Speed Range:	750 - 3,200 RPM
Time:	1 - 60 Minutes or Continuous
Vortexing Orbit:	1.0mm Diameter
Maximum Stirring Capacity:	1 Microplate with Standard SBS Footprint
Dimensions (D x W x H):	6.5 x 7.0 x 5.5 Inches
Weight:	6.1lbs

Part Number	Description	List Price	SALE Price
CLS-3064-001	Microplate Genie, 120V	525.00	325.00*