

LIMITED OFFER!

Only ONE at this price!
Show Model with Platform



Only 11" wide.

\$1800 for Both



**MAGic Clamp Universal Platform
and Rubber Mat Included!**

(Clamps and flasks are not included.)

INCU-SHAKER™ Mini Incubator with Rubber Mat and MAGic Clamp™ Universal Platform

- Turned "on" but never used
- Tradeshow display unit ♦
- 1 year warranty

This Incu-Shaker™ Mini combines digitally controlled temperature and circular (orbital) mixing, to optimize the cultivation of a wide variety of cell types. It's also the most compact shaking incubator in its class. Despite the remarkably small footprint, a large 11.5 x 9.5in workspace accepts a variety of platforms for use with flasks, tubes and other common culture vessels.

The standard platform is a non-slip, rubber-coated surface, ideal for tissue culture flasks, Petri dishes and staining trays.

This limited offer also includes a MAGic Clamp™ Universal platform that provides a suitable surface for use with MAGic clamps, the easiest way to instantly change between flask clamps of different sizes. Each clamp utilizes a powerful (patent pending) magnetic base to securely attach the clamp to the shaking platform without the use of screws, bolts or tools. Clamps are not included and need to be purchased separately. **See the second page of this flyer for a listing of MAGic Clamps.**

Temperature, shaking speed and run-time are displayed on the large LED control panel located on the front of the unit. The integral microprocessor includes a constant monitoring system that verifies and maintains accuracy through the duration of the program. Sophisticated over-temperature and over-speed controls ensure long life, safety and sample integrity.

Part Number	Description	Reg. Prices	Clearance Prices
CLS-2543-001	Mini Incubator, Rubber Mat, 120V	\$2,650.00 *	\$1,800.00 *
CLS-2543-003	Universal Platform	209.00 *	Included!
		\$2,859.00*	\$1,800.00 *

Save \$1059.00

**Use promo code:
CG-14061**



Where chemistry meets biology.

Call or email to place your order: 1-800-843-1794 • customer-service@cglifesciences.com

* Not Subject to Discount. Trademarks are the property of their respective owners.

♦Tradeshow display equipment may show minor wear or surface imperfections due to handling.

MAGic Clamps and Accessories for MAGic Clamp Platforms

Part Number	Description	Max. # of Items		Price
		Incu-Shaker Mini		
CLS-2543-050	Erlenmeyer Clamps for MAGic Platform, 50mL	20		25.00 *
CLS-2543-125	Erlenmeyer Clamps for MAGic Platform, 125mL	12		26.00 *
CLS-2543-250	Erlenmeyer Clamps for MAGic Platform, 250mL	8		28.00 *
CLS-2543-500	Erlenmeyer Clamps for MAGic Platform, 500mL	5		29.50 *
CLS-2543-1L	Erlenmeyer Clamps for MAGic Platform, 1000mL	2		48.00 *
CLS-2543-016	Microplate Clamps for MAGic Platform	4		68.50 *
CLS-2543-005	Conical Tube Rack, 12 x 50mL and 21 x 15mL	1		129.00 *
CLS-2543-017	Conical Tube Rack, 32 x 15mL	2		129.00 *
CLS-2543-018	Conical Tube Rack, 15 x 50mL	2		129.00 *
CLS-2543-019	Conical Bottle Rack, 4 x 500mL or 600mL	1		149.00 *



* Not Subject to Discount.

Call or email to place your order: 1-800-843-1794 • customer-service@cglifesciences.com