

Chemglass

MICROWAVE

Vials



2-5mL 10-20mL 0.5-2mL 0.2-0.5mL

NEW Packaging!

- Designed for convenience and product protection.
- Ten vials per removable tray.
- Stir bars included.



Updated Blue Septa!



MICROWAVE VIAL PACKS

Microwave vials are manufactured from heavy wall Type I, Class A borosilicate glass to withstand pressures up to 300 PSI (20 bars). The 20mm aluminum seals are supplied with 0.125" thick blue PTFE / white silicone septa installed. For crimper and decapper tools see CG-4930.

Part No.	Description	OD x Length (mm)	Qty. Included in Pack			Price
			Vials	Seals	Stir Bars	
CG-4920-01	2-5mL Vials	16 x 83	100	100	100	390.00 †
CG-4920-02	10-20mL Vials	28 x 83	100	100	12	1,061.10 †
CG-4920-03	0.5-2mL Tapered Vials	16 x 83	100	100	100	490.00 †
CG-4920-04	0.2-0.5mL Vials	8 x 83	100	100	100	825.30 †

Replacement Parts		Case Qty.	Price
CG-4920-10	20mm Silver Aluminum Seals with 0.125" thick Blue PTFE/White Silicone Septa	200	190.55 †
CG-4920-15	Stir Bars, PTFE, 3mm Dia. x 10mm Long	50	88.30 †
CG-4920-16	Stir Bars, PTFE, 6mm Dia. x 15mm Long, Egg-Shaped	12	22.70 †
CG-4920-17	Stir Vanes, PTFE, 5.25mm Wide x 11mm High	50	134.20 †
CG-4920-18	Stir Bars, PTFE, 2mm Dia. x 5mm Long	50	107.75 †

† Subject to Reusable Glass / Plastic Discount. CG-15035IHPDF0515



Phone 1-800-843-1794
1-856-696-0014

Fax 1-800-922-4361
1-856-696-9102

customer-service@cglifesciences.com
www.cglifesciences.com