



Vacuum Adapters

CG-1049-A ADAPTER, LUER, VACUUM

Adapter has an outer luer joint that provides a convenient way to connect solid phase extraction (SPE) cartridges, tubes or syringes to a standard taper inner joint. The adapter provides for easier sample concentration and clean-up while reducing the chances of cross-contamination. Hose connection is 10mm OD at the largest serration.

Part Number	Description	Price
CG-1049-A-14	Adapter, Luer Outer to 14/20 Inner	56.10
CG-1049-A-24	Adapter, Luer Outer to 24/40 Inner	58.31

CG-1049-B ADAPTER, LUER, VACUUM, STOPCOCK

Adapter has an outer luer joint that provides a convenient way to connect solid phase extraction (SPE) cartridges, tubes or syringes to a standard taper inner joint. The adapter provides for easier sample concentration and clean-up while reducing the chances of cross-contamination. The addition of a 2mm Teflon stopcock between the outer luer joint and standard taper inner permits control of the vacuum source. Hose connection is 10mm OD at the largest serration.

Part Number	Description	Price
CG-1049-B-14	Adapter, Luer Outer to 14/20 Inner	83.45
CG-1049-B-24	Adapter, Luer Outer to 24/40 Inner	85.66

CG-1049-F ADAPTER, FEMALE LUER LOCK, VACUUM

Adapter has a 1/4" -28 thread that accepts a female luer lock adapter which provides a convenient way to connect solid phase extraction (SPE) cartridges, tubes, or syringes with a male luer lock fitting to a standard male taper inner joint. The adapter provides for easier sample concentration and clean-up while reducing the chances of cross-contamination. Hose connection is 10mm OD at the largest serration. Supplied complete with 1 each CG-1055-20; 1/4"-28 Female to Female luer lock adapter.

Part Number	Description	Price
CG-1049-F-14	Adapter, Luer Lock Female to 14/20 Inner	82.25
CG-1049-F-24	Adapter, Luer Lock Female to 24/40 Inner	85.20

Replacement Components:

Part Number	Description	Price
CG-1055-20	1/4"-28 Female to Female Luer Lock Adapter	12.85*

CG-1049-G ADAPTER, FEMALE LUER LOCK, VACUUM, STOPCOCK

Adapter has a 1/4" -28 thread that accepts a female luer lock adapter which provides a convenient way to connect solid phase extraction (SPE) cartridges, tubes, or syringes with a male luer lock fitting to a standard male taper inner joint. The adapter provides for easier sample concentration and clean-up while reducing the chances of cross-contamination. The addition of a 2mm Teflon stopcock between the luer joint and standard taper inner permits control of the vacuum source. Hose connection is 10mm OD at the largest serration. Supplied complete with 1 each CG-1055-20; 1/4"-28 Female to Female luer lock adapter.

Part Number	Description	Price
CG-1049-G-14	Adapter, Luer Lock Female to 14/20 Inner	110.25
CG-1049-G-24	Adapter, Luer Lock Female to 24/40 Inner	113.00

Replacement Components:

Part Number	Description	Price
CG-1055-20	1/4"-28 Female to Female Luer Lock Adapter	12.85*

*not subject to discount



Vacuum Adapters



CG-1052-A ADAPTER, VACUUM FILTRATION

Adapter used for reduced pressure filtration with plain stem buchner funnels. Top has a HEAVY DUTY tooled flare to accept the listed size pluro stopper. Hose connection is 10mm OD at the largest serration. ADAPTER IS NOT SUPPLIED WITH RUBBER PLURO STOPPER.

<u>Part Number</u>	<u>Joint Size</u>	<u>Flange ID (mm)</u>	<u>Pluro Stopper</u>	<u>Overall Height (mm)</u>	<u>Accepts up to Buchner Funnel Capacity (mL)</u>	<u>Price</u>
CG-1052-A-14	14/20	23	#2	51	60	46.98
CG-1052-A-24	24/25	32	#3	65	150	53.89



PLURO STOPPERS:

<u>Part Number</u>	<u>Adapter Size</u>	<u>Price</u>
CG-1401-11	#2	4.65*
CG-1401-12	#3	5.13*

*not subject to discount



3800 North Mill Road • Vineland NJ 08360
 U.S. and Canada phone: 1-800-843-1794 • Fax: 1-800-922-4361
 customer-service@cglifesciences.com • www.cglifesciences.com

CHEMGLASS LIFE SCIENCES IS A WOMAN-OWNED BUSINESS