

2010
PRICING

Disposable smartSpatulas®**


Used for sample handling and manipulation in labs across the country.

These versatile spatulas can be used to scoop, scrape or stir. Made of virgin polypropylene, smartSpatulas are an economic alternative to reusable spatulas.

Available in 3 convenient sizes: Macro-Green 310 mm, Standard-Blue and Standard-Opaque 210 mm, and Antistatic-Micro 140 mm.

smartSpatulas are autoclavable for sterile applications. Single-use smartSpatulas help to reduce cross contamination, are resistant to dilute acids and bases, and the unique hollow core minimizes thermal transfer.

smartSpatulas comply with FDA regulations for food-contact applications, are shatter resistant for liquid nitrogen frozen samples and can also be used as a stir rod or policeman.



“if you use disposable weigh boats... you **Should** use disposable smartSpatulas!”

INEXPENSIVE



SHATTER-RESISTANT



ANTISTATIC MICRO



- INEXPENSIVE spatulas that SCOOP, SCRAPE & STIR.
- AUTOCLAVABLE for sterile applications.
- ANTISTATIC MICRO for small volumes.
- Unique hollow cores MINIMIZE THERMAL TRANSFER.
- Excellent for use as STIR ROD OR POLICEMAN.
- Resistant to dilute acids & bases
- Comply with FDA regulations for food-contact applications. Made of virgin PP.
- SHATTER RESISTANT for liquid nitrogen frozen samples.

©2005, 2006, 2007 Investigen, Inc. Patent Pending. smartSpatula is a trademark of Investigen, Inc. All rights reserved.

PART NUMBER	SPATULA LENGTH	TYPE/COLOR	CASE SIZE	PRICE/PKG. <small>*NOT SUBJECT TO DISCOUNT</small>
CG-1985-P-01	210mm	Standard, Blue	300/Pack	\$26.67*
CG-1985-P-02	210mm	Standard, Blue	7500/Case	\$666.67*
CG-1985-P-03	210mm	Standard, Opaque	300/Pack	\$26.67*
CG-1985-P-04	210mm	Standard, Opaque	7500/Case	\$666.67*
CG-1985-P-05	140mm	Antistatic, Micro, Clear	300/Pack	\$22.33*
CG-1985-P-06	140mm	Antistatic, Micro, Clear	7500/Case	\$492.14*
CG-1985-P-07	310mm	Macro, Green	150/Pack	\$28.99*
CG-1985-P-08	310mm	Macro, Green	3600/Case	\$724.64*



3800 North Mill Road • Vineland NJ 08360 • U.S. and Canada phone: 1-800-843-1794
Fax: 1-800-922-4361 • customer-service@cglifesciences.com • www.cglifesciences.com