

# CHEMGLASS

## PIE-BLOCK™

Please Call for  
Current Pricing.

CG-1991-P REACTION, PIE-BLOCK™,  
SAFETY HOLDER

- Interchangeable Pie Wedges Permit Various Configurations On One Hot Plate
- Pie Wedges Available for:
  - ▲ Vials - 1 dram to 40mL
  - ▲ 10mL to 100mL RBF
  - ▲ Microwave Vials
- Inserts Are Color Coded For Quick And Easy Identification
- Holes Are Spaced For Optimum Stirring Of All Vials
- Holder Fits All 135mm Diameter Hot Plate Stirrers

New Sizes  
Available



Safety Pie Plate Holder  
Keeps Up To Four Pie Wedges  
Securely Locked Together  
On Hotplate

Aluminum Construction  
Provides Excellent  
Heat Transfer  
While Not Interfering  
With The Magnetic Field

Ordering Information:  
1-800-843-1794

Pie-Block™ is a Trademark of Chemglass, Inc.

# CHEMGLASS

3800 North Mill Road • Vineland NJ 08360  
U.S. and Canada phone: 1-800-843-1794 • Fax: 1-800-922-4361  
customer-service@chemglass.com • www.chemglass.com