

CHEMGLASS

PIE-BLOCK™

CG-1991-P REACTION, PIE-BLOCK™,
SAFETY HOLDER

Please Call for
Current Pricing

New Sizes
Available

- Interchangeable Pie Wedges Permit Various Configurations On One Hot Plate
- Pie Wedges Available for:
 - ▲ Vials - 1 dram to 40mL
 - ▲ 10mL to 100mL RBF
 - ▲ Microwave Vials
- Inserts Are Color Coded For Quick And Easy Identification
- Holes Are Spaced For Optimum Stirring Of All Vials
- Holder Fits All 135mm Diameter Hot Plate Stirrers

Safety Pie Plate Holder
Keeps Up To Four Pie Wedges
Securely Locked Together
On Hotplate

Aluminum Construction
Provides Excellent
Heat Transfer
While Not Interfering
With The Magnetic Field

Ordering Information:
1-800-843-1794

Pie-Block™ is a Trademark of Chemglass, Inc.

CHEMGLASS

3800 North Mill Road • Vineland NJ 08360
U.S. and Canada phone: 1-800-843-1794 • Fax: 1-800-922-4361
customer-service@chemglass.com • www.chemglass.com

CG-1991-P REACTION, PIE-BLOCK™, SAFETY HOLDER

CHEMGLASS
PIE-BLOCK

NEW Reaction PIE-BLOCK™ is designed to fit directly on the CG-1990 or CG-1990-A circular top hot plate stirrers (135mm Diameter). Versatile design of the safety pie plate holder allows up to four pie wedges on one hot plate and securely locks them together. All wedges are anodized aluminum and color coded for quick and easy identification.

A drilled thermowell is provided for insertion of the CG-1995-25 Opti-Chem® Electronic Contact Thermometer or CG-1995-02 PT-100 Temperature Sensor.

Complete PIE-BLOCK™ system supplied with one safety pie plate holder and one each of the six wedges. (Codes -10 through -15)

*Custom configurations are available.
Please call for pricing.*



PART

NUMBER

DESCRIPTION

PRICE

CG-1991-P-01	PIE-BLOCK™ System, Supplied With 1 Safety Pie Plate Holder With 1 each of the Six Listed Pie Wedges (Code -10 thru -15)	316.02*
--------------	---	---------

REPLACEMENT PIE PLATE SAFETY HOLDER:

PART

NUMBER

DESCRIPTION

PRICE

CG-1991-P-05	Reaction Pie Plate Safety Holder Only, Green Anodized	58.92*
--------------	---	--------

VIAL PIE WEDGES:

PART NUMBER

VIAL TYPE

VIAL SIZE (mm)

ANODIZED COLOR

HOLE DEPTH (mm)

OF VIALS

PRICE EACH

CG-1991-P-10	20mL (Scintillation)	28 x 60	Black	24	4	44.99*
CG-1991-P-11	30mL (Scintillation)	28 x 73	Green	30	4	49.28*
CG-1991-P-12	40mL (Scintillation)	28 x 98	Orange	43	4	53.14*
CG-1991-P-13	1 Dram (4mL)	15 x 45	Gold	19	9	55.17*
CG-1991-P-14	2 Dram (8mL)	17 x 60	Blue	25	8	53.29*
CG-1991-P-15	4 Dram (16mL)	21 x 70	Red	32	4	49.28*
CG-1991-P-19	2-5mL (Microwave)	-	Torquoise	-	8	55.89*
CG-1991-P-20	10-20mL (Microwave)	-	Berry	-	4	53.82*
CG-1991-P-23	.5-2mL (Microwave)	-	Purple	-	8	67.28*

FLASK PIE WEDGES:

PART NUMBER

FLASK SIZE (mL)

ANODIZED COLOR

OF FLASKS

PRICE EACH

CG-1991-P-30	10mL	Gold	1	\$66.00*
CG-1991-P-32	25mL	Blue	1	\$67.50*
CG-1991-P-34	50mL	Red	1	\$69.00*
CG-1991-P-36	100mL	Black	1	\$70.50*

* Not subject to discount

CHEMGLASS

3800 North Mill Road • Vineland NJ 08360
U.S. and Canada phone: 1-800-843-1794 • Fax: 1-800-922-4361
customer-service@chemglass.com • www.chemglass.com