

CG-1994

CHEMGLASS OPTIMAG+, MAGNETIC HOT PLATE STIRRER WITH SAFETY CONTROL

NEW!



Our new OPTIMAG+ hot plate stirrer with safety control features. This economically priced unit has many of the same features as our more expensive units.

FEATURES:

- Set Safety Temperature Limit Displayed Digitally
- Exact Temperature And Speed Setting Via Digital Display, Even When Switched Off
- Hot Top Indicator >> Hot Surface Warning To Prevent Burns!
- Indicator Lights For Both Heating And Stirring Modes
- Enclosed Housing Prevents Liquids And Vapors From Damaging Internal Components
- Highly Polished Aluminium Heating Plate For Optimum Heat Transfer

SPECIFICATIONS:

H2O Stirring Quantity MAX (L).....	20
Plate Diameter (mm):.....	135
Stirring Speed (RPM):.....	50 to 1,500
Temperature Range (°C):.....	RT-310
Overall Dimensions, W x D x H (mm):.....	165 x 275 x 85
Electrical.....	115V, 60 Hz
Weight (lbs):.....	5.5

PART

NUMBER	DESCRIPTION	PRICE
CG-1994-10	Hot Plate Stirrer with PT1000 Sensor	640.00*
CG-1994-50	OPTIMAG+ Magnetic Hot Plate Stirrer Package	685.00*

COMPLETE PACKAGE TO INCLUDE 1 EACH:

CG-1990-A-15	PT 1000 Temperature Sensor
CG-1995-05	Boss Head
CG-1995-06	S.S. Vertical Support Rod
CG-1995-22	S.S. Holding Rod

*Not Subject to Discount



LIT-120104-01