

LABeler™ Handheld Laboratory Printer



- **Lightweight, handheld operation**
- **Variety of label tapes for many tube sizes**
- **Built in cutter - no need for scissors**
- **Barcoding** 
- **“Ultra-Low” tape for liquid nitrogen storage**

Identifying and organizing increasingly large numbers of samples has become an important objective for today's busy research laboratory. The process of writing on sample tubes & containers can be messy and unreliable. It can also be difficult to find the time to access a computer and design a layout for expensive ink jet labels. The new LABeler™ Lab Printer provides the convenience of a lab marking pen with the benefits of a large desktop printer. Simply type in the desired information using the 50+ button key pad and press the print key...just press the cut lever and your label is instantly ready for application.


With over 12 available sizes, colors & material combinations, the LABeler is ideal for label identification of any laboratory vessel for nearly any laboratory application. Standard label tapes can be autoclaved, boiled or even frozen to temperatures as low as -50°C. Special Ultra-Low label tapes are also available for storage in -80°C environments or for liquid nitrogen applications.



Available with 9V power supply
(also runs off 6 AA batteries)



Replaceable tape cassettes
(26ft / 792cm length)



www.cglifesciences.com

Ordering Information

Cat. No.	Description	Price
CG-4034-10	LABeler™ Handheld lab printer, with 9V power supply (U.S. Plug)*	\$149.95*
CG-4034-12	LABeler™ Handheld lab printer, without power supply	\$109.95*
CG-4034-14	Power supply (9V) for LABeler™ Lab Printer (U.S. Plug)*	\$41.50*
CG-4034-16	26' Cassette of 6mm lab tape, white with black print	\$19.50*
CG-4034-18	26' Cassette of 6mm lab tape, clear with black print	\$24.50*
CG-4034-20	26' Cassette of 12mm lab tape, white with black print	\$22.00*
CG-4034-22	26' Cassette of 12mm lab tape, clear with black print	\$27.50*
CG-4034-24	26' Cassette of 12mm lab tape, yellow with black print	\$22.00*
CG-4034-26	26' Cassette of 12mm lab tape, red with white print	\$22.00*
CG-4034-28	26' Cassette of 12mm lab tape, blue with white print	\$22.00*
CG-4034-30	26' Cassette of 12mm lab tape, green with white print	\$22.00*
CG-4034-32	26' Cassette of 24mm lab tape, white with black print	\$28.00*
CG-4034-34	26' Cassette of 24mm lab tape, clear with black print	\$35.00*
CG-4034-36	20' Cassette of 12mm ULT (ultra-low temp) lab tape, white with black print	\$25.00*
CG-4034-38	20' Cassette of 24mm ULT (ultra-low temp) lab tape, white with black print	\$35.00*

* Life Science Equipment Product Line

Toll Free: 1-800-843-1794

email: customer-service@cglifesciences.com

www.cglifesciences.com