



CG-4901-150 PRE-PACKED HPLC VIALS, 15MM X 45MM

Chemglass has made it easier for you by providing convenient pre-packed vials. Every pack of 15mm Dia. x 45mm long threaded vials comes complete with the corresponding stir bar already inside each vial. Threaded closures are also supplied in separate bag. Vials are packaged in convenient LAB-PACS of 144 pieces.

COMPLETE PACKAGES:

PART

<u>NUMBER</u>	<u>DESCRIPTION</u>	<u>PRICE</u>
CG-4901-150	15mm x 45mm Clear Vial, 4mL, 13-425 Threaded Top, Supplied with 8mm x 1.5mm Dia. Stir Bar & Screw Cap with PTFE/Silicone Face Septa, Pack of 144	\$389.00†



CG-4901-125 PRE-PACKED HPLC VIALS, 12MM X 32MM

Chemglass has made it easier for you by providing convenient pre-packed vials. Every pack of 12mm Dia. x 32mm long threaded vials comes complete with the corresponding stir bar already inside each vial. Threaded closures are also supplied in separate bag. The components are packaged in clear trays with clear snap on lids. The trays fit easily into bench drawers for easy access and the packaging keeps products visible & particle free.

COMPLETE PACKAGES:

PART

<u>NUMBER</u>	<u>DESCRIPTION</u>	<u>PRICE</u>
CG-4901-125	12mm x 32mm Clear Vial, 2mL, Graduated, 10-425 Threaded Top, Supplied with 7mm x 2.0mm Dia. Stir Bar & Screw Cap with PTFE/Red Rubber Face Septa, Pack of 100	\$268.50†



CG-4901-120 PRE-PACKED HPLC VIALS, 12MM X 32MM

Chemglass has made it easier for you by providing convenient pre-packed vials. Every pack of 12mm Dia. x 32mm long crimp top vials comes complete with the corresponding stir bar already inside each vial. Crimp top closures & Septa are also supplied in separate bag. The components are packaged in clear trays with clear snap on lids. The trays fit easily into bench drawers for easy access and the packaging keeps products visible & particle free.

COMPLETE PACKAGES:

PART

<u>NUMBER</u>	<u>DESCRIPTION</u>	<u>PRICE</u>
CG-4901-120	12mm x 32mm Clear Vial, 2mL, Graduated, Silver Aluminum Crimp Top, Supplied with 7mm x 2.0mm Dia. Stir Bar & Screw Cap with PTFE/Red Rubber Face Septa, Pack of 100	\$251.10†

† subject to consumable discount