

NEW

**STAINLESS
STEEL**

LABORATORY CLAMPS & CONNECTORS



CG-9205-S CLAMP, THREE-PRONG EXTENSION, STAINLESS STEEL

These durable, multi-purpose clamps are made entirely of stainless steel with an electro-polished finish offering excellent chemical resistance and overall strength. The versatile 3-prong design secures various lab apparatus and the dual adjustments offer a wide range of motion. Supplied with vinyl and fiberglass sleeve prong covers for high temperature applications.

Part Number	Size	Min. to Max Grip Size	Arm Length	Arm Dia.	Overall Length	Price
CG-9205-S-04	Small	0 to 48mm	102mm	10mm	172mm	60.90*
CG-9205-S-05	Med.	0 to 69mm	127mm	13mm	229mm	70.50*
CG-9205-S-06	Large	0 to 102mm	127mm	13mm	273mm	78.00*

*not subject to discount

CG-9252-S-01 CLAMP HOLDER, LARGE, STAINLESS STEEL

This heavy duty holder is made of stainless steel with an electro-polished finish with greater grip strength and the ability to hold flat panels or shelves up to 3/4" (19mm) thick. Large thumb screws allow easy adjusting.



<u>Part Number</u>	<u>Min. to Max. Grip Size</u>	<u>Price</u>
CG-9252-S-01	0 to 19mm (3/4")	41.70*

CG-9262-S-01 HOOK CONNECTOR, STAINLESS STEEL

This easy-to-use connector has a stainless steel electro-polished finish and allows one-handed assembly of two components with one adjustment screw. Maximum diameter of 1/2" (13mm) for each component.



<u>Part Number</u>	<u>Price</u>
CG-9262-S-01	23.50*



<u>Part Number</u>	<u>Price</u>
CG-9264-S-01	24.60*



<u>Part Number</u>	<u>Price</u>
CG-9274-S-01	37.70*

CG-9264-S-01 S-CONNECTOR, STAINLESS STEEL

This clamp has a stainless steel electro-polished finish and features ergonomic contouring for easy adjustment with angled set screws for maximum grip against rods or clamps.

CG-9274-S-01 LAB FRAME FOOT, STAINLESS STEEL

This strong, durable lab-frame foot has a stainless steel electro-polished finish which is designed to permanently mount to bench tops, floors or walls. Spilled fluids roll right off the smooth body contour, minimizing contamination. Accommodates 1/2" (13mm) diameter rod. Supplied with three screws.