

NEW



CG-1761 VIAL FUNNEL, TEFLON

CG-1761-20

Vial Funnel Kit, Supplied With A Polypropylene Storage Box And One Of The Eight Funnel Sizes

\$99.⁰⁰*

Replacement Components

Part Number	Funnel OD (in)	G.P.I. Thread Size	Price <small>*NOT SUBJECT TO DISCOUNT</small>
CG-1761-01	1.0	8-425	\$11.18*
CG-1761-02	1.0	9-425, 10-425	\$11.18*
CG-1761-03	1.0	13-425	\$11.39*
CG-1761-04	1.4	15-425	\$13.46*
CG-1761-05	1.4	18-400	\$13.46*
CG-1761-06	1.9	20-400	\$15.94*
CG-1761-07	1.9	22-400	\$15.94*
CG-1761-08	1.9	24-400	\$15.94*

Cross Reference Chart

Vial Type & Capacity	Vial OD X Length (mm)	Thread Size
2mL Standard Opening	12 x 32	8-425
2mL Big Mouth Vial	12 x 32	10-425
1 Dram (4mL)	15 x 45	13-425
1.5mL CD™	13 x 44	13-425
5mL CD™	15 x 51	13-425
4mL WISP™	15 x 45	13-425
2 Dram (8mL)	17 x 60	15-425
5mL CD™	17 x 51	15-425
4 Dram (16mL)	21 x 70	18-400
10mL CD™	21 x 60	18-400
20mL CD™	28 x 67	22-400
20mL Scintillation	28 x 60	24-400
30mL Scintillation	28 x 73	24-400
40mL Scintillation	28 x 98	24-400
40mL CD™	28 x 103	24-414
60mL CD™	28 x 143	24-414

Please Call
1-800-843-1794
For Ordering
Info. & Current
Pricing

*not subject to discount



CG-4900

VIAL, SAMPLE, CLEAR, TEFLON LINED CAP

Clear sample vials packaged in convenient LAB-PACS of 144 pieces. Made from Type I borosilicate glass. Full height partitions minimize breakage. Solid Teflon lined caps. Open top caps with septa are available under part number CG-4910.

Part Number	Vial Size OD x Height (mm)	Vial Capacity (mL)	Screw Thread Size	Price/Pkg.
CG-4900-01	15 x 45	4	13-245	\$105.57
CG-4900-02	17 x 60	8	15-425	\$111.78
CG-4900-03	21 x 70	16	18-400	\$121.10
CG-4900-04	28 x 57	20	24-400	\$178.02
CG-4900-05	28 x 70	30	24-400	\$183.20
CG-4900-06	28 x 95	40	24-400	\$198.72

CG-4905

VIAL, SAMPLE, AMBER, TEFLON LINED CAP

Amber sample vials packaged in convenient LAB-PACS of 144 pieces. Made from Type I borosilicate glass. Full height partitions minimize breakage. Solid Teflon lined caps. Open top caps with septa are available under part number CG-4910.



Part Number	Vial Size OD x Height (mm)	Vial Capacity (mL)	Screw Thread Size	Price/Pkg.
CG-4905-01	15 x 45	4	13-245	\$134.55
CG-4905-02	17 x 60	8	15-425	\$144.90
CG-4905-03	21 x 70	16	18-400	\$156.29
CG-4905-04	28 x 57	20	24-400	\$232.88
CG-4905-05	28 x 70	30	24-400	\$241.16
CG-4905-06	28 x 95	40	24-400	\$259.79

◀ ◀ Ordering Information 1-800-843-1794 ▶ ▶