

**30%  
COST  
SAVINGS**



# ROTARY EVAPORATOR BUNDLES

**30%  
COST  
SAVINGS**



**FEATURING** the IKA RV10 Digital Rotovap with motorized lift, digital RPM display and a heating bath with digital temperature control

» CG-1335-R-10



» 21 Lpm, 1.5 torr with silencer



## CGS-1335-R80

### Vacuum Pump Package

- Compact and quiet
- Suitable for rotary evaporation of high boiling point solvents
- Service interval of 10,000-15,000 hours (5 years of typical use)

**\$5,720.00**

» CG-1335-R-10



» Vacuum Controller and Solenoid Valve



» 21 Lpm, 1.5 torr with silencer



## CGS-1335-R82

### Vacuum Pump/Controller Package - all features of the above, plus

- Electronic two-point control and integrated vent valve
- Up to ten multi-step user programs
- Easy-to-use controller with large display and graphic interface

**\$7,895.00**

» CG-1335-R-10



» 27 Lpm, 1.5 torr



## CGS-1335-R84

### Digital Self-Adjusting Vacuum Pump System Package - all features of above, plus

- Just press "Start" automated control for automatic boiling point determination with continuous optimization, for up to 30% faster evaporation
- Integrated solvent recovery
- Virtually silent with reduced energy consumption, and greatly extended service intervals

**\$8,495.00**