

## CGS-1990-140

### IKAMAG® RCT Basic Magnetic Hot Plate Stirrer KIT

The **NEW** CG-1990-14 RCT Basic IKAMAG® magnetic hot plate stirrers incorporate an adjustable safety circuit and are indisputably at the top of their class in worldwide technology. The front panel design allows easy control and reading of both the speed and temperature features. The heating plate is 135mm diameter and made from highly polished Silumin for optimum heat transfer. The IKAMAG® hot plate stirrers are protected against liquid and vapor penetration by an enclosed assembly (IP 42) which guarantees long service life. The overhanging corner protects the instrument's plugs and electrical connections from spills.

✓ **COMPLETE KIT INCLUDES:** Hot Plate Stirrer, CG-1995-25 ETS D6 Sensor, Vertical Support Rod, Boss Head, S.S. Holding Rod.

✓ **FREE WITH EACH KIT:** 1ea CG-1991-01 REACTION BLOCK, 20mL Scintillation Vial: Black anodized aluminum reaction blocks offer excellent heat transfer while not interfering with the magnetic stirrer.

\*Vials must be purchased Separately

**\$920<sup>00</sup>**  
SALE PRICE



✓ **FREE**  
**REACTION**  
**BLOCK**  
**WITH KIT!**