


# High Vacuum Valves Precision Bore



## CGK-560 & 534 >> 1-ARM KONTES® STYLE VALVE

This Kontes® style hi-vac valve has a single sidearm seal at a 90° angle to the barrel. This inexpensive configuration allows for the lower portion of the barrel to be sealed directly to the apparatus being constructed, making for more compact, stronger assemblies. Designed with a precision external thread for fine control. The CGK-534 has an exposed o-ring tip, while the CGK-560 has a teflon tip. Fully interchangeable with Kontes® style valves.



Part Number	Description	Price
CGK-534-01	0-4mm 1-Arm, Exposed O-Ring Tip	47.96
CGK-534-02	0-8mm 1-Arm, Exposed O-Ring Tip	76.94
CGK-534-03	0-12mm 1-Arm, Exposed O-Ring Tip	81.00
CGK-560-01	0-4mm 1-Arm, Teflon Tip	44.90
CGK-560-02	0-8mm 1-Arm, Teflon Tip	76.94
CGK-560-03	0-12mm 1-Arm, Teflon Tip	81.00

### Replacement Components

Part Number	Description	Price
CGK-502-01	0-4mm 1-Arm, Glass Only	30.61
CGK-502-02	0-8mm 1-Arm, Glass Only	35.97
CGK-502-03	0-12mm 1-Arm, Glass Only	46.85
CGK-535-01	0-4mm Plug Only, Exposed O-Ring Tip	18.37
CGK-535-02	0-8mm Plug Only, Exposed O-Ring Tip	41.84
CGK-535-03	0-12mm Plug Only, Exposed O-Ring Tip	44.00
CGK-561-01	0-4mm Plug Only, Teflon Tip	15.31
CGK-561-02	0-8mm Plug Only, Teflon Tip	41.84
CGK-561-03	0-12mm Plug Only, Teflon Tip	44.00

## CGK-551 & 525 >> 90° KONTES® STYLE VALVE

Hi-vac valve with two sidearms sealed to the barrel at a 90° angle to one another. Designed with a precision external thread for fine control. The CGK-525 has an exposed o-ring tip, while the CGK-551 has a teflon tip. Fully interchangeable with Kontes® style valves.



Part Number	Description	Price
CGK-525-01	0-4mm 90°, Exposed O-Ring Tip	56.12
CGK-525-02	0-8mm 90°, Exposed O-Ring Tip	87.86
CGK-525-03	0-12mm 90°, Exposed O-Ring Tip	92.50
CGK-551-01	0-4mm 90°, Teflon Tip	53.06
CGK-551-02	0-8mm 90°, Teflon Tip	87.76
CGK-551-03	0-12mm 90°, Teflon Tip	92.50

### Replacement Components

Part Number	Description	Price
CGK-503-01	0-4mm 90°, Glass Only	36.73
CGK-503-02	0-8mm 90°, Glass Only	54.08
CGK-503-03	0-12mm 90°, Glass Only	56.70
CGK-535-01	0-4mm Plug Only, Exposed O-Ring Tip	18.37
CGK-535-02	0-8mm Plug Only, Exposed O-Ring Tip	41.84
CGK-535-03	0-12mm Plug Only, Exposed O-Ring Tip	44.00
CGK-561-01	0-4mm Plug Only, Teflon Tip	15.31
CGK-561-02	0-8mm Plug Only, Teflon Tip	41.84
CGK-561-03	0-12mm Plug Only, Teflon Tip	44.00

## CGK-529 & 555 >> 0-4MM, IN-LINE KONTES® STYLE VALVE

Hi-vac valve with bent arms to minimize hold up of fluids that would occur in conventional valves. Designed with a precision external thread for fine control. The CGK-529 has an exposed o-ring tip, while the CGK-555 has a teflon tip. Fully interchangeable with Kontes® style valves.



Part Number	Description	Price
CGK-529-01	0-4mm In-Line, Exposed O-Ring Tip	76.53
CGK-555-01	0-4mm In-Line, Teflon Tip	73.47

### Replacement Components

Part Number	Description	Price
CGK-505-01	0-4mm In-Line, Glass Only	57.14
CGK-535-01	0-4mm Plug Only, Exposed O-Ring Tip	18.37
CGK-561-01	0-4mm Plug Only, Teflon Tip	15.31

## CGK-527 & 553 >> S-TYPE KONTES® STYLE VALVE

Hi-vac valve with two sidearms sealed to the barrel at a 180° angle to one another and on the same plane. The CGK-527 has an exposed o-ring tip, while the CGK-553 has a teflon tip. Fully interchangeable with Kontes® style valves.



Part Number	Description	Price
CGK-527-01	0-4mm 180°, Exposed O-Ring Tip	78.57
CGK-527-02	0-8mm 180°, Exposed O-Ring Tip	92.04
CGK-527-03	0-12mm 180°, Exposed O-Ring Tip	120.00
CGK-553-01	0-4mm 180°, Teflon Tip	75.51
CGK-553-02	0-8mm 180°, Teflon Tip	92.04
CGK-553-03	0-12mm 180°, Teflon Tip	120.00

### Replacement Components

Part Number	Description	Price
CGK-504-01	0-4mm 180°, Glass Only	61.22
CGK-504-02	0-8mm 180°, Glass Only	72.45
CGK-504-03	0-12mm 180°, Glass Only	94.00
CGK-535-01	0-4mm Plug Only, Exposed O-Ring Tip	18.37
CGK-535-02	0-8mm Plug Only, Exposed O-Ring Tip	41.84
CGK-535-03	0-12mm Plug Only, Exposed O-Ring Tip	44.00
CGK-561-01	0-4mm Plug Only, Teflon Tip	15.31
CGK-561-02	0-8mm Plug Only, Teflon Tip	41.84
CGK-561-03	0-12mm Plug Only, Teflon Tip	44.00

## CGK-531 & 557 >> 0-4MM, 3-WAY KONTES® STYLE VALVE

Hi-vac valve with three way configuration for use in construction of manifolds, adapters, etc. The CGK-531 has an exposed o-ring tip, while the CGK-557 has a teflon tip. Fully interchangeable with Kontes® style valves.



Part Number	Description	Price
CGK-531-01	0-4mm 3-Way, Exposed O-Ring Tip	127.55
CGK-557-01	0-4mm 3-Way, Teflon Tip	124.49

### Replacement Components

Part Number	Description	Price
CGK-506-01	0-4mm 3-Way, Glass Only	110.20
CGK-535-01	0-4mm Plug Only, Exposed O-Ring Tip	18.37
CGK-561-01	0-4mm Plug Only, Teflon Tip	15.31

### Replacement Components

Standard for all valves shown here



**CGK-535-01**  
0-4mm Plug Only, Exposed O-Ring Tip  
**CGK-535-02**  
0-8mm Plug Only, Exposed O-Ring Tip  
**CGK-535-03**  
0-12mm Plug Only, Exposed O-Ring Tip



**CGK-561-01**  
0-4mm Plug Only, Teflon Tip  
**CGK-561-02**  
0-8mm Plug Only, Teflon Tip  
**CGK-561-03**  
0-12mm Plug Only, Teflon Tip