

Replacement Glassware for Automated Reactors

Compatible with Mettler® Reactors†




CG-MR-50V1

• **REACTION VESSEL, FLAT BOTTOM, 50mL, METTLER STYLE**

74.5mm OAH, 45mm OD Flat Flange.

Equivalent to Mettler part number 51161252.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-50V1	FB Reaction Vessel, 50mL, Mettler Style	164.00*



CG-MR-50V2

• **REACTION VESSEL, ROUND BOTTOM, 50mL, METTLER STYLE**

74.5mm OAH, 45mm OD Flat Flange.

Equivalent to Mettler part number 51161255.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-50V2	RB Reaction Vessel, 50mL, Mettler Style	163.00*



CG-MR-250V

• **REACTION VESSEL, ROUND BOTTOM, 250mL, METTLER STYLE**

140mm OAH, 79mm OD Body, 98mm OD O-Ring Flange.

Equivalent to Mettler part number 51161059.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-250V	RB Reaction Vessel, 250mL, Mettler Style	360.00*

• **MICRO COIL CONDENSER, METTLER STYLE**

140mm OAL, 14/20 top outer joint, 14/20 lower inner joint, GL-14 thread inlet/outlet, Supplied complete with caps with apertures and hose connections.

Equivalent to Mettler part number 51161029.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-250C1	Micro Coil Condenser, Mettler Style	125.00*

• **MICRO COIL CONDENSER, JACKETED, METTLER STYLE**

160mm OAL, 14/20 top outer joint, 14/20 lower inner joint, GL-14 thread inlet/outlet, Supplied complete with caps with apertures and hose connections.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-250C2	Jacketed Micro Coil Condenser, Mettler Style	233.00*



CG-MR-250C1



CG-MR-250C2

*not subject to discount

†Mettler is a registered trademark of Mettler Toledo



CG-MR-50L1

LID, REACTION VESSEL, 4-NK, 50mL, METTLER STYLE

(1) 14/20 outer SN, (1) GL-14 SN, (1) 10/18 outer SN, (1) 10/10 outer SN, 45mm OD Flange.

Equivalent to Mettler part number 51161032.


<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-50L1	Lid, Reaction Vessel, 4-Nk, 50mL, Mettler Style	305.00*

LID, REACTION VESSEL, 6-NK, 50mL, METTLER STYLE

(1) 14/20 outer SN with 10/18 outer, (2) GL-14 SN's, (1) 10/18 outer SN, (1) 7/12 outer or (1) 7/15 outer SN and 45mm OD Flange.

Equivalent to Mettler part number 51161100 and 51160817.

<u>Part Number</u>	<u>Description</u>	<u>Overall Height (A)</u>	<u>Price</u>
CG-MR-50L2	Lid, Reaction Vessel, 6-Nk	73mm	310.00*
CG-MR-50L3	Lid, Reaction Vessel, 6-Nk	63mm	315.00*



CG-MR-250L 2

LID, SHORT, REACTION VESSEL, 6-NK, 250mL, METTLER STYLE

14/20 outer CN, (1) 14/20 outer SN's, (2) 19/22 outer SN, and (2) GL-14 SN's, 98mm OD Flange with notch.

Equivalent to Mettler part number 51161188.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-250L1	Lid, Reaction Vessel, 6-Nk, 250mL, Mettler Style	370.00*

LID, REACTION VESSEL, 6-NK, 250mL, METTLER STYLE

14/20 outer CN, (3) 14/20 outer SN's, and (2) 19/22 outer SN's, 98mm OD Flange without notch.

Equivalent to Mettler part number 51161072.

<u>Part Number</u>	<u>Description</u>	<u>Price</u>
CG-MR-250L2	Lid, Reaction Vessel, 6-Nk, 250mL, Mettler Style	370.00*

*not subject to discount



CG-MR-250L1



CG-MR-250L2



3800 North Mill Road
 Vineland NJ 08360
 U.S. and Canada phone:
 1-800-843-1794
 Fax: 1-800-922-4361
 customer-service@cglifesciences.com
 www.cglifesciences.com