



REACTION VIALS

CG-4904 REACTION VIAL, 20mL, 28 x 57mm

Screw cap vials are designed for heated sample reactions.

Product Highlights:

- ▶ **Packaged in a clear PETG tray complete with a slide or lift off dust cover.** Vials are separated by a SBS partition to eliminate glass to glass contact. Caps are supplied complete in the package, not attached. Entire package is recyclable.
- ▶ **Vials are made from 33 expansion borosilicate glass** which conforms to ASTM E-438 Type 1-Class A, USP Type 1 and ISO 3585. The high resistance to thermal shock along with its maximum chemical resistance makes this the glass of choice for reaction vials along with most all laboratory glassware.
- ▶ **Polypropylene caps have a septa bonded to the cap,** using the Chemglass SURE-LINK™ process. This inert process bonds the .125" Teflon Faced Silicone septa directly to the cap without adhesives, making the caps autoclavable and able to withstand multiple injections.
- ▶ **All vials are manufactured with the larger 24-400 thread,** providing easier access for sampling or automation work.
- ▶ **Available with a choice of caps:** green open top caps or white solid top caps with SURE-LINK™ Teflon faced septa. The white caps are ideal for storage purposes.



Part Number	Description	Price
CG-4904-01	Reaction Vial, 20mL, Complete with Green Open Top Caps, TFE Septa	\$65.80▲
CG-4904-03	Reaction Vial, 20mL, Complete with White Solid Caps, TFE Septa	\$71.80▲

▲ Subject to the Vials/Containers Discount

STORAGE/REACTION VIALS

CG-4908 STORAGE/REACTION VIAL, 20mL, 28 x 57mm

For storage applications, these vials are packaged in a SBS coated white tray complete with SBS partition. Caps are supplied installed on the vials. Available with a choice of WHITE solid top or GREEN open top caps, both having a TFE lined silicone SURE-LINK™ septa.

Part Number	Description	Price
CG-4908-01	Reaction Vial, 20mL, Complete with Green Open Top Caps, TFE Septa	\$58.80▲
CG-4908-03	Reaction Vial, 20mL, Complete with White Solid Caps, TFE Septa	\$64.80▲

CG-4909 STORAGE/REACTION VIAL, 40mL, 28 x 95mm

Similar to CG-4908, but 40mL instead of 20mL. These vials are packaged in a SBS coated white tray complete with SBS partition. Caps are supplied installed on the vials. Available with a choice of open top or solid top caps, both having a TFE lined silicone SURE-LINK™ septa.

Part Number	Description	Price
CG-4909-05	Reaction Vial, 40mL, Complete with Green Open Top Caps, TFE Septa	\$89.00▲
CG-4909-07	Reaction Vial, 40mL, Complete with White Solid Caps, TFE Septa	\$95.00▲

CG-4911 SURE-LINK™ REPLACEMENT CAPS

Polypropylene caps have a septa bonded to the cap, using the Chemglass SURE-LINK™ process. This inert process bonds the .125" Teflon Faced Silicone septa directly to the cap without adhesives, making the caps autoclavable and able to withstand multiple injections. The white caps are ideal for storage purposes.

Part Number	Description	Price
CG-4911-20	Cap, Green, Open Top, SURE-LINK™ Septa, PKG of 100	\$50.00▲
CG-4911-21	Cap, White, Solid Top, SURE-LINK™ Septa, PKG of 100	\$56.00▲

▲ Subject to the Vials/Containers Discount



CG-4908

Packaging printed with alpha numeric coordinates for easy identification.



CG-4909



CG-4911-20



CG-4911-21