



REACTION VIALS 1 AND 2 DRAM

CG-4904 REACTION VIAL, 1 AND 2 DRAM

Screw cap vials are designed for heated sample reactions.

Product Highlights:

- ▶ **Packaged in a clear PETG tray complete with a slide or lift off dust cover.** Vials are separated by a SBS partition to eliminate glass to glass contact. Caps are supplied complete in the package, not attached. Entire package is recyclable.
- ▶ **Vials are made from 33 expansion borosilicate glass** which conforms to ASTM E-438 Type 1-Class A, USP Type 1 and ISO 3585. The high resistance to thermal shock along with its maximum chemical resistance makes this the glass of choice for reaction vials along with most all laboratory glassware.
- ▶ **Polypropylene caps have a septa bonded to the cap,** using the Chemglass Life Sciences SURE-LINK™ process. This inert process bonds the Teflon Faced Silicone septa directly to the cap without adhesives, making the caps autoclavable and able to withstand multiple injections.
- ▶ **Available with a choice of caps:** green open top with SURE-LINK™ or green solid top caps.



REACTION VIAL, 1 AND 2 DRAM

OPEN TOP:

Part Number	Description	Capacity	Thread Finish	Pkg. Qty.	Price
CG-4904-05	Reaction Vial, 1 Dram, 15 x 45mm, Complete with Green Open Top Polypropylene Cap, TFE/Silicone Septa	4mL	13-425	100	\$56.00▲
CG-4904-08	Reaction Vial, 2 Dram, 17 x 60mm, Complete with Green Open Top Polypropylene Cap, TFE/Silicone Septa	8mL	15-425	100	\$61.00▲



CG-4904-05
CG-4904-08

SOLID TOP:

Part Number	Description	Capacity	Thread Finish	Pkg. Qty.	Price
CG-4904-06	Reaction Vial, 1 Dram, 15 x 45mm, Complete with Green Solid Top Cap, TFE Septa	4mL	13-425	100	\$49.00▲
CG-4904-09	Reaction Vial, 2 Dram, 17 x 60mm, Complete with Green Solid Top Cap, TFE Septa	8mL	15-425	100	\$56.00▲



CG-4904-06
CG-4904-09

REPLACEMENT CAPS ONLY:

Part Number	Description	Pkg. Qty.	Price
CG-4911-10	13-425 Solid Top Caps	100	\$36.00▲
CG-4911-12	15-425 Solid Top Caps	100	\$37.00▲
CG-4910-15	13-425 Open Top SURE-LINK™ Caps	100	\$39.00▲
CG-4910-16	15-425 Open Top SURE-LINK™ Caps	100	\$40.00▲



CG-4911-10
CG-4911-12



CG-4910-15
CG-4910-16

▲ Subject to the Vials/Containers Discount

* SURE-LINK™ is a trademark of Chemglass Life Sciences



3800 N. Mill Road, Vineland, NJ 08360 | US and Canada (P) 1-800-843-1794 (F) 1-800-922-4361
OUTSIDE US AND CANADA (P) 856-696-0014 (F) 856-696-7366 | customer-service@cglifesciences.com
www.cglifesciences.com