

Tissue Grinders

These tissue grinders are engineered for safety and ease of use. The patented closed system design provides the ability to grind the specimen in a sealed container helping to prevent any risk of contact with the sample. The narrow pestle eliminates the problem of sample overflow common with other tissue grinders. The molded-in abrasive surface does not leave any sediment to obstruct the sample during examination. Tissue grinders come fully assembled and gamma sterilized in an easy-to-open package. An extra cap and patient identification label are provided for storage, transport and traceability.

Part Number	Description	Case Weight	Case Qty.	Case Price
CLS-5015	15mL, Conical	3 lbs.	10	92.47 [◇]
CLS-5035	35mL, Free Standing	4 lbs.	10	92.47 [◇]
CLS-5050	50mL, Conical	4 lbs.	10	92.47 [◇]

[◇] Subject to LS consumable discount



ULTRA Tissue Grinders

The pestle tips of the new Ultra Tissue Grinders have been vitrified using a proprietary silica casting process. The result is a glass-like abrasive surface which provides a more aggressive grinding action to enable efficient processing of difficult and fibrous samples. Ultra Tissue Grinders incorporate the same patented closed system as the standard grinders. Available in 15mL and 50mL sizes, they are supplied gamma sterilized in an easy-to-open package with an extra cap for storage and transport plus a patient identification label for traceability.

Part Number	Description	Case Weight	Case Qty.	Case Price
CLS-5115	15mL, Conical	3 lbs.	10	114.81 [◇]
CLS-5150	50 mL, Conical	4 lbs.	10	114.81 [◇]

[◇] Subject to LS consumable discount



ULTRA Vitrified Tip For Grinding

Difficult Samples:

- Muscle
- Connecting Tissue

Standard Molded-In Abrasive Tip