



CG-1992 OptiTherm® Quick Setup Guide

*U.S. Patent # 7,494,267



Setup for Digital Thermoregulator:

Option 1

Chemglass hotplate with optional Digital Thermoregulator inserted into OptiTherm block thermowell. (Thermowells are designed for plain 1/8" diameter thermocouples.)

CAUTION
Maximum Temperature
+180°C



Option 2

Chemglass hotplate with optional Digital Thermoregulator inserted into flask.



Setup for PT100 Sensor:

Option 1:

Chemglass hotplate with PT100 Sensor inserted into OptiTherm block thermowell. (Thermowells are designed for plain 1/8" diameter thermocouples.)



CAUTION
Maximum Temperature
+180°C

Option 2:

Chemglass hotplate with PT100 Sensor inserted into flask.

