

DATALOGGER

dual function temperature recorder and data logger

Easy to use handheld device takes and records temperature readings via an attached probe. A large touch screen, color display indicates temperature, maximum and minimum readings and will produce a graph of the data.



In 'Data Logging' mode (Not CFR21 compliant.)

- The unit can be set up and left to take temperature readings at predetermined intervals for a desired length of time

In 'Measurement' mode, the unit can be used in two ways

- Predefined Checklist - used as an electronic checklist to do **spot temperature readings** on a predefined list of items
- Spot Measurement - used as a thermometer to capture **one-off temperature readings** displayed on screen

This data can then be downloaded to PC via a USB connection and opened with Microsoft Excel® or the software available on our website. The unit is compatible with any J, K and T type thermocouple probes.

Note: Probes and cables must be purchased separately.

Large 2.8" (71 mm) color touch screen display



Additional probe styles and sizes available on our website

Part No.	Description	Price
CG-3497-100	Digital Temperature Monitor with Data Logging	\$225.00 †
Extension Cords		
CG-3499-01	10ft Ext. Cord, Type "J" (Black), Standard to Mini Plug	\$44.50 *
CG-3499-02	10ft Ext. Cord, Type "T" (Blue), Standard to Mini Plug	\$40.20 *
CG-3499-03	10ft Ext. Cord, Type "K" (Yellow), Standard to Mini Plug	\$36.10 *

†Reusable Glass / Plastics *Not Subject To Discount.



1-800-843-1794
www.cglifesciences.com
tech@cglifesciences.com

Excel is a Registered Trademark of Microsoft
 Printed in USA. CG-17058BC1.5M617



See other side for details on unique water-saving Findenser